

A Combinatorial Hydrogel Platform to Probe Stem Cell Chondrogenesis in 3D

Sebastián L. Vega, Kwang Hoon Song, Chao Wang, Lin Han, Jason A. Burdick

Monday October 24, 2016, 2:05PM



The Stem Cell Microenvironment



Challenges in deciphering cell-ECM interactions:


- Tedious to evaluate many conditions
- Unable to decouple mechanical and biochemical cues

Need for novel materials and high-throughput techniques to probe cell-ECM interactions



Techniques to probe cell-biomaterial interactions

ECM Mimetic Arrays




ECM Mimetic Arrays

- Complex DOE
- Expensive equipment not readily available

1st Generation Combinatorial Hydrogels

- Mostly explored in 2D contexts
- 2D does not represent 3D environments

Combinatorial Hydrogels




Dolatshahi-Pirouz et al., *Sci Rep*, 2014;
Rape et al., *Nat Commun*, 2015.

Combinatorial hydrogels to screen 3D cell-ECM interactions


Design Criteria:

- 3D cell encapsulation
- Cytocompatible crosslinking
- Ability to vary biochemical and mechanical signals
- Rapid single cell imaging possible

Thiol-ene reactions between norbornene and thiols




Scheme for fabricating combinatorial hydrogels




Gramlich et al., *Biomaterials*, 2013.

RGD gradients to probe effects on chondrogenesis

Chondrogenic Differentiation




Experimental Design



Vinatier et al., Trends in Biotech, 2009.



Confocal microscopy of 5x5 mm hydrogels

Develop whole-hydrogel imaging protocol




Develop Sox9 and aggrecan image analysis protocols

Sox9 (+1 day)





Aggrecan (+7 days)






Combinatorial hydrogel characterization

Hydrogel characterization
• RGD gradient




- Tile-scan confocal imaging
- Measure horizontal fluorescence (x1 → x10)
- 5 horizontal regions (n=3)
- Irradiate gels, measure fluorescence of known RGD concentrations



Combinatorial hydrogel characterization


Hydrogel characterization

- Mechanics



7.1 ± 1 kPa


Hertz indentation model



10x10 array, 15 measurements/position

$$F = \frac{4}{3} \cdot \frac{E_s}{1 - \mu^2} \cdot \sqrt{r} \cdot \sigma^{\frac{3}{2}}$$

E_s = hydrogel modulus
 μ = Poisson's ratio
 r = probe radius



Elastic moduli (5 – 9 kPa)

7.8	8.4	8.6	8.6	8.0	8.9	7.9	7.9	8.5	6.6
8.2	8.9	8.6	7.1	8.6	8.9	8.3	8.9	8.6	8.6
8.1	8.4	8.4	7.9	8.8	8.5	8.2	6.6	6.2	5.7
7.8	7.6	7.6	7.9	7.6	8.1	6.5	5.8	5.8	6.0
6.9	7.9	7.9	7.4	8.0	6.7	6.3	5.8	6.7	5.6
7.8	7.9	8.5	6.7	6.0	6.0	5.7	6.8	5.3	5.6
7.8	7.2	5.9	7.1	7.6	6.5	6.5	5.8	5.7	5.6
5.9	7.2	7.2	7.5	7.2	7.6	6.3	6.3	5.8	5.6
7.2	6.9	7.5	7.8	6.5	6.3	6.4	6.2	5.6	5.5
7.4	7.9	6.7	7.2	7.5	7.0	6.2	6.0	5.2	5.6

Variation from mean (< 30%)




0.09	0.17	0.21	0.20	0.12	0.25	0.11	0.11	0.19	0.07
0.14	0.25	0.21	0.01	0.20	0.25	0.16	0.25	0.21	0.21
0.14	0.17	0.18	0.10	0.23	0.19	0.15	0.08	0.13	0.20
0.10	0.06	0.07	0.11	0.06	0.14	0.10	0.19	0.19	0.15
0.03	0.11	0.11	0.04	0.12	0.06	0.12	0.19	0.07	0.22
0.09	0.10	0.19	0.06	0.15	0.15	0.21	0.05	0.26	0.21
0.10	0.00	0.17	0.00	0.06	0.09	0.08	0.18	0.21	0.21
0.17	0.01	0.01	0.05	0.01	0.06	0.11	0.12	0.19	0.22
0.01	0.03	0.05	0.10	0.09	0.12	0.10	0.14	0.21	0.23
0.04	0.11	0.06	0.01	0.04	0.02	0.13	0.16	0.27	0.22

Assessment of MSCs encapsulated in combinatorial hydrogels

Cellular Studies



- Divide gel into 10 vertical regions (x1 → x10)
- Divide each region into 10 squares
- For each square:
 - **Viability:** count live/dead cells
 - **Sox9:** nuclear fluorescence
 - **Aggrecan:** secreted volume





Viability

- High MSC viability for at least 7 days in combinatorial hydrogels
- Rapid imaging & image analysis developed (>300 cells analyzed per position)

Assessment of MSCs encapsulated in combinatorial hydrogels



Nuclear Sox9 fluorescence (+1d)

•No RGD (x10) → highest nuclear Sox9




Aggrecan production (+7d)

•Best: x8 (~1mM RGD), x10 (no RGD)




Summary

Developed combinatorial hydrogel system




- High viability
- Wide range of biochemical cues (0 → 5mM)
- No significant change in elastic moduli


Rapid confocal imaging and image analysis



Significant differences in chondrogenic markers



Can add additional biochemical cues



Acknowledgments and questions

Polymeric Biomaterials Lab

Dr. Jason A. Burdick

Dr. Kwang Hoon Song

Sara Trujillo-Muñoz

John Bricker

Michelle Kwon

Dr. Claudia Loebel

Renata Lang Sala

Dr. Chris Highley

Dr. Yi-Cheun Yeh

Dr. Adrienne Rosales

Dr. Dawn Pedrotty

Dr. Su Chin Heo

Dr. Steven Caliarì

Dr. Kalil Abdullah

Dr. Chris Rodell

Andrew Rouff

Minna Chen

Josh Mealy

Leo Wang

Liliang Ouyang

James Howard

Jonathan Galarraga

Mu-Huan Lee



NanoBiomechanics Lab

Dr. Lin Han

Chao Wang

